Aim-T_{conv} HMP

For human Pan T cell
Activation and Expansion
FOR RESEARCH USE ONLY

Version 1.1.2

Other materials required:

- ☐ Cryopreserved or freshly isolated T Cells, or CD3⁺ T Cells
- ☐ RPMI 1640 supplemented with 2 mM L-glutamine, or equivalent culture medium.
- ☐ Heat inactivated Fetal Bovine Serum (FBS)
- ☐ Penicillin/streptomycin
- ☐ Recombinant human IL-2 (rIL-2)
- ☐ 2-Mercaptoethanol
- ☐ Cell culture vessels
- ☐ Humidified CO₂ incubator or bioreactor

Part A) Activation and expansion:

1. **Medium preparation:** RPMI + 10% FBS + 1% Pen/Strep + 20~50 μ M 2-Mercaptoethanol + 30 to 100 U/mL rIL-2, or other compatible culture media.
2. **HMP preparation:** Resuspend Aim-T_{conv} hydrogel microparticles (HMP) in the vial by vortexing.
3. **HMP seeding:** Calculate desired HMP seeding volume per well. Freshly prepare adequate HMP dilutions in culture medium for easy and accurate pipetting. Aliquot diluted HMP in to each well.
 - Recommended T cell seeding density : 2 x 10⁵ to 3 x 10⁵ cell per cm²
 - Recommended HMP to PBMC ratio: 2:1.
 - Recommended HMP to T cell ratio: 1:1.

Example T cell seeding plan				HMP vol. to add	
Plate	Area (cm ²)	T Cell/ well	T Cell/ cm ²	1x HMP ratio	2x HMP ratio
96-well	0.32	8 x 10 ⁴	2.5 x 10 ⁵	2 μ L	4 μ L
48-well	0.95	2 x 10 ⁵	2.1 x 10 ⁵	5 μ L	10 μ L
24-well	1.9	4 x 10 ⁵	2.1 x 10 ⁵	10 μ L	20 μ L
6-well	9.5	2 x 10 ⁶	2.1 x 10 ⁵	50 μ L	100 μ L
T25	25	6 x 10 ⁶	2.4 x 10 ⁵	150 μ L	300 μ L
T75	75	2 x 10 ⁷	2.7 x 10 ⁵	500 μ L	1000 μ L

4. **T Cell seeding:** Aliquot resuspended T cells. Gently mix T cells and HMP by pipetting up/down 3 to 5 times, ensure HMP and cells are evenly distributed under microscope.
5. **Co-culture:** Incubate in a humidified 5% CO₂ incubator at 37°C. Monitor T cell morphology and confluency every other day. **DO NOT** disturb HMP-cell interaction in the first 4 days.
6. **HMP cleanup:** HMP will self-degrade by hydrolysis after around 14 days under culture environment. Alternatively, HMP can be rapidly enzymatically degraded by adding the included HMP digesting buffer directly into the culture medium at 1X final concentration. After incubating at 37°C for 1-2 hours, the HMPs will be completely degraded
7. **Restimulation:** Restimulation is required for expansion longer than 10 days. T cell growth usually slows down after day 10 after initial activation. Recommended HMP to cell ratio for restimulation is 0.5:1 to 1:1.

Website

allegrowbiotech.com

LinkedIn

allegrow-biotech

Email

welovecells@allegrowbiotech.com

Component Descriptions:

- A suspension of hydrogel microparticles (HMP) made of chemically crosslinked dextran.
- A vial of HMP digesting buffer (10X, 1.0mL), containing enzymes to digest HMP.
- HMPs sized 6-15 μ m were coated with phospholipid bilayer, with membrane docked T cell activating signal panels.
- Supplied at 4 x 10⁷ / mL, in 0.5 mL PBS per vial.
- **Stability and Storage:**
- Shipped with ice, keep product refrigerated (2 to 8°C).
- Stable at 4°C for 6 months. Contents are sterile in unopened tube.
- Do not subject product to freezing, high temperature (greater than 40°C).

Things to note:

When seeding T cells with HMP

Ensure T cells and HMP are evenly distributed to maximize interaction.

When you see T cell cluster.

HMP and T cells tend to aggregate in the well center over time. Gently shake the culture plate to redistribute HMP and T cells. **Avoid unnecessary pipette mixing in the first 4 days**, disturbing the HMP-cell clustering will cause suboptimal cell growth.

When to add new medium or split?

Monitor T cell growth periodically by performing T cell sampling and counting periodically. Supplement fresh medium or pass the cells to new culture vessels when

- Color turns orange yellow (acidic, ~ pH 6.5).
- T Cell grows to over 3 x 10⁶ cells /mL.

When to restimulate?

Day 8 – 10 after previous stimulation.


Part B) T Cell phenotype characterization

1. Typical T cell morphology:
 - Resting T cells: smaller in size, round shaped.
 - Activated T cells: are larger in size, irregular.
2. Characterize T cells with FACS accordingly:

	Marker	Type
Lineage	CD3+	Pan T
	CD4+	CD4 T
	CD8+	CD8 T
	CD45RA+ CD62L+	T _{naive}
	CD45RA- CD62L+	T _{central memory}
	CD45RA+ CD62L-	T _{effector memory}
	CD45RA- CD62L-	T _{effector}
	CD45RO+	T _{cm} & T _{em}
Activation	CD25	Activated T
Exhaustion	PD-1	Exhausted T


REPRESENTATIVE DATA

Human Pan T cells co-cultured with Aim-T_{conv} at different time points


Representative bright-field micrographs showing the activation and expansion of Pan T cells at various time points. Scale bar = 50 μm. Pan T cells were activated and expanded using Aim-Tconv HMP at 1:1 HMP: T cell ratio. Red signals indicate the HMPs.


1 T cell proliferation profile


2 T cell viability post-activation

3 T cell activation profiles

4 T cell memory subtypes

1. T Cell growth (n=4 T cells from healthy human donors) after activation by Aim-Tconv HMP and other popular T cell activation reagents. Expansion fold = viable cells at time t / initial number of T cells.
2. T Cell viability (n=4 T cells from healthy human donors) after activation by Aim-Tconv and other popular T cell activation reagents.
3. The expression of early stage marker (CD69) and late stage marker (CD25) indicating T cell activation. Unstimulated T cell was used as gating.
4. T cell memory subtypes after activation by Aim-Tconv and other popular T cell activation reagents.

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