



## Bst DNA Polymerase

### User's Instruction

#### Description

Bst DNA Polymerase is a homologous protein of Bst DNA Polymerase, large fragment, derived from the thermophilic bacterium *Bacillus stearothermophilus* (Bst). When compared to the large fragment of Bst DNA Polymerase, it exhibits stronger 5'→3' DNA polymerase activity, enhanced strand displacement capability, tolerance to dUTP, salt resistance, and resistance to non-ionic detergents. Bst DNA Polymerase lacks 5'→3' and 3'→5' exonuclease activity. It can be used in various applications, including loop-mediated isothermal amplification (LAMP), crossing priming amplification (CPA), rolling-circle amplification (RCA), and isothermal amplification reactions based on rolling-circle amplification.

The isothermal amplification temperature mediated by Bst DNA Polymerase generally falls between 50-68°C, typically at 65°C. The optimal temperature depends on the primers and the products being amplified and may require experimental optimization.

Compared to the **Bst 2.0 DNA Polymerase** from similar companies, this product exhibits **similar enzyme activity, comparable high dUTP tolerance, similar high tolerance to non-ionic detergents, and equivalent salt resistance.**

#### Kit Contents

	40KU
1. Bst DNA Polymerase (40U/μl)	1 ml
2. 10 × Isothermal Buffer	3 ml
3. 100 mM MgSO <sub>4</sub>	2 ml

#### Protocol

- Set up isothermal amplification reaction as the following table:

Component	Volume	
Nuclease-free Water	(15.6-x)μl	-
	1	



10 × Isothermal Buffer	2.5μl	1X
MgSO <sub>4</sub> (100mM)	1.5μl	6mM (8mM total)
dNTP (25mM each)	1.4μl	1.4mM each
FIP/BIP Primers (25X, 40μM)	1μl	1.6μM
F3/B3 Primers (25X, 5μM)	1μl	0.2μM
Loop F/B Primers (25X, 10μM)	1μl	0.4μM
Template	xμl	> 10 copies or more
Bst DNA Polymerase (40U/μl)	1μl	1600U/ml
Total volume	25μl	-

#### Note

- After the completion of the reaction system setup, add an appropriate amount (1μl) of high-concentration SYBR Green I to the reaction tube cap for every 25μl of the reaction system. After the isothermal amplification reaction is finished, centrifuge at 8000×g for 1 minute. A positive reaction is indicated by the fluorescence turning green, while a negative reaction is indicated by the reaction system remaining colorless or brown. Alternatively, an indicator may not be added; a positive reaction can be recognized by the reaction solution becoming visibly turbid after the reaction program, whereas a negative reaction is indicated by the reaction solution remaining clear.
- For reaction optimization, adjust the Mg<sup>2+</sup> concentration (4-10mM), enzyme amount (0.04-0.32U/μl), or change the reaction temperature (50-68°C).
- If analysis is conducted through agarose gel electrophoresis or any method requiring the opening of the LAMP reaction container, establish separate areas and equipment for auxiliary analysis to prevent contamination.
- Due to the rapid nature of the reaction, to ensure reproducibility of the experiment, it is recommended to add the template DNA at the final stage.
- It is strongly recommended to include a negative control without the template to ensure the specificity of the amplification.
- To prevent contamination during reagent preparation, it is essential to operate within a clean bench.
- The preparation of reagents and template DNA should preferably be conducted in different areas from the analysis of PCR products by techniques such as gel electrophoresis to avoid contamination.

#### 2. Thermo Condition:

- a) 65°C for 60 min



b) 85°C for 20 min (inactivation)

## Storage

Store at -20°C. Avoid multiple freeze-thaw cycles.

**Only for research and not intended for treatment of humans or animals**