

**NEVER  
STOP**

**FUJIFILM**  
Value from Innovation

## FUJIFILM COVID-19 Ag Test

### IVDD Technical Documentation File Results in Performance Evaluation

Version 3



## **Materials about performance characteristics**

### **for FUJIFILM COVID-19 Ag Test**

#### **1. Specimen type**

##### **Appropriate specimen type(s)**

Nasopharyngeal swabs

##### **Specimen Stability**

No data

##### **Freeze/thaw cycle**

No data

##### **Preparation of sample**

Collected specimen should be prepared as a sample in accordance with the proper method.

## 2. Performance ( TC ST-0202-SIL-COV-001 )

### (i) Accuracy

#### 1) PROTOCOL

In accordance with Procedure, each specimen was measured for 3 reagent lots.

#### 2) CRITERIA

- Testing an in-house positive control sample<sup>1)</sup> results in a positive result.
- Testing an in-house negative control sample<sup>2)</sup> results in a negative result.

Note 1) A sample made from SARS-CoV-2 recombinant protein diluted to 25 pg/mL in Extraction Reagent Solution R1

Note 2) Extraction Reagent Solution R1.

#### 3) RESULT

Specimen	Lot	Result
in-house negative control	900001	Negative
	900002	Negative
	900003	Negative
in-house positive control	900001	Positive
	900002	Positive
	900003	Positive

#### 4) CONCLUSION

Satisfied criteria and acceptable.

## (ii) Precision

### 1) PROTOCOL

In accordance with Procedure, each specimen was measured 3 times at the same time for 3 reagent lots.

### 2) CRITERIA

- Testing three in-house positive control samples<sup>1)</sup> at the same time yields all positive results.
- Testing three in-house negative control samples<sup>2)</sup> at the same time yields all negative results.

Note 1) A sample made from SARS-CoV-2 recombinant protein diluted to 25 pg/mL in Extraction Reagent Solution R1

Note 2) Extraction Reagent Solution R1.

### 3) RESULT

Specimen	Lot	Result
in-house positive control	900001	Positive Positive Positive
	900002	Positive Positive Positive
	900003	Positive Positive Positive
in-house negative control	900001	Negative Negative Negative
	900002	Negative Negative Negative
	900003	Negative Negative Negative

### 4) CONCLUSION

Satisfied criteria and acceptable.

**(iii) Minimum detection sensitivity (Limit of Detection) (Aras No.:001-021101)**

**1) PROTOCOL**

A sample collected from a patient suspected of being infected with SARS-CoV-2 at the Yokohama City University Hospital, and a nasopharyngeal swab storage sample diagnosed as SARS-CoV-2 negative by RT-PCR was used as a pool sample. Inactivated SARS-CoV-2 (JPN / TY / WK-521 strain) was added to this pool sample, and prepared dilution. In accordance with Procedure, each dilution was measured.

**2) RESULT**

Lot.900002

Specimen	Lot	Concentration [TCID <sub>50</sub> /mL]	Result
Inactivated SARS-CoV-2 (JPN / TY /WK-521 strain)	900002	2.3 × 10 <sup>2</sup>	Positive Positive Positive
		1.1 × 10 <sup>2</sup>	Positive Positive Positive
		5.6 × 10 <sup>1</sup>	Negative Positive Negative
		2.8 × 10 <sup>1</sup>	Negative Negative Negative

**3) CONCLUSION**

The limit of detection is concentration equivalent to inactivated SARS-CoV-2 (JPN / TY / WK-521 strain) 1.1×10<sup>2</sup> TCID<sub>50</sub>/mL.

### 3. Correlation (Aras No.:001-021140)

SARS-CoV-2 infection positive / negative of each patient was judged by reference method;

Positive: RT-PCR for nasopharyngeal swabs testing positive.

Negative: RT-PCR for nasopharyngeal swabs testing negative.

#### 1) PROTOCOL

In accordance with Procedure, each nasopharyngeal swabs specimen was measured and compared to the determination result of the reference method.

#### 2) RESULT

Lot.900002 and Lot. 210120

PCR-Positive

		RT-PCR Ct value			
		< 26	< 28	< 30	< 35
FUJIFILM COVID-19 Ag Test	Positive	100	107	110	110
	Negative	0	3	6	11
	Total	100	110	116	121
	Sensitivity (PPA)	100.0%	97.3%	94.8%	90.9%
	95% CI	96.4-100.0%	92.2-99.4%	89.1-98.1%	84.3-95.4%

PCR-Negative

		RT-PCR
		Negative
FUJIFILM COVID-19 Ag Test	Positive	0
	Negative	104
	Total	104
	Specificity (NPA)	100.0%
	95% CI	96.5-100.0%

Overall results of clinical performance

		RT-PCR		
		Positive	Negative	Total
FUJIFILM	Positive	110	0	110
COVID-19 Ag	Negative	11	104	115
Test	Total	121	104	225

Sensitivity (PPA):90.9%(110/121, 95% CI:84.3-95.4%)

Specificity (NPA):100.0%(104/104, 95% CI:96.5-100.0%)

**EXPLANATION OF TERMS:**

Sensitivity or Positive Percent Agreement (PPA) =

True Positives / (True Positives + False Negatives)

Specificity or Negative Percent Agreement (NPA) =

True Negatives / (True Negatives + False Positives)

**3) CONCLUSION**

Good correlation result was obtained with the reference method.

#### **4. Interfering substances and medications (TC ST-0206-SIL-COV-001)**

##### **1) PROTOCOL**

Prepare a sample by adding each interfering substance to in-house negative control, in-house positive control.

In accordance with Procedure, each sample was measured.

##### **2) RESULT**

Lot.900003

At the concentrations listed below, the following substances will not invalidate the results given by this product.

- Hemoglobin from hemolysis (450 mg/dL)
- Blood (2%)



## 5. Cross reactivity (TC ST-0205-SIL-COV-001)

### 1) PROTOCOL

Prepare a sample by adding each recombinant coronavirus antigens and inactivated influenza virus to in-house negative control, in-house positive control.

In accordance with Procedure, each sample was measured.

### 2) RESULT

Lot.900003

No cross reactivity was observed with the following virus.

<Recombinant coronavirus antigens>

- MERS-CoV
- HCoV-229E,
- HCoV-OC43
- HCoV-NL63
- HCoV-HKU1
- SARS-CoV

<Inactivated influenza virus>

- Influenzavirus H1N1
- Influenzavirus H3N2
- Influenzavirus B

# FUJIFILM



**FUJIFILM Corporation**

26-30, NISHIAZABU 2-CHOME, MINATO-KU, TOKYO 106-8620, JAPAN



**FUJIFILM Europe GmbH**  
Heesenstrasse 31, 40549 Duesseldorf, Germany

EU Importer:  
**FUJIFILM Europe B.V.**  
Oudenstaart 1, 5047 TK Tilburg, The Netherlands