

# **Chemiluminescence (CL) Application Note**

## Phosphatidylcholine Hydroperoxide of Increased in Type 2 Diabetic Patients

#### **Aims**

This study was designed to clarify serum phosphatidylcholine hydroperoxide levels related to blood glucose control in Type 2 diabetic patients.

### Method

HPLC column: SIL-NH2

Mobile phase:

2-propanol-methanol-water (135:45:20, v/v/v)

CL reagent:


10 mg/L Cytochrome c and 1.0 mg/L luminol in borate

buffer (pH 10.0)

Standard:

Photo-oxidized L-  $\alpha$  -phosphatidylcholine,  $\beta$  -oleoyl-  $\gamma$  -palmitoyl (C18:1,[cis] -9/C16:0, SIGMA)

Detector: CLA-2100, CLA-FL2



#### Result

The data showed that (1) serum phosphatidylcholine hydroperoxide is increased in Type 2 diabetes patients, (2) the level strongly correlates with HbA<sub>1c</sub> rather than the fasting blood glucose level.

Oxidative stress is markedly observed in Type 2 diabetes patients with an increase of serum PCOOH. Such oxidative stress might trigger the chronic complications in diabetes mellitus.

### Reference

Nagashima, T., Oikawa, S. and Yamada, R. : Diabetes Research and Clinical Practice, in press




Fig. 2 Phosphatidylcholine hydroperoxide (PCOOH) levels in serum-PCOOH when normal volunteers and diabetic patients were subgrouped into 6 groups according to  $HbA_{1c}$  levels: (1) normal volunteers, (2) less than 6, (3) 6-6.4, (4) 6.5-6.9, (5) 7.0-7.4, and (6) over than 7.5 %.



URL: http://www.tei-c.com

No. clbe2014\_004v4

Contact us sales@toi-s.con

TOHOKU ELECTRONIC INDUSTRIAL CO.,LTD

Contact us... sales@tei-c.com