



VIOLET & BLUE RAY PROTOCOL GUIDE





WELCOME

Welcome to the user guide for your new Violet or Blue Ray device! This device is a type of electrotherapy that has been used for over a century to treat a variety of health conditions. Violet or Blue ray devices work by producing a high-frequency electrical current that is applied to the skin via a glass electrode. The current and combined argon and neon gases create a Violet or Blue-colored glow in the electrode, which is where the device gets its name.

In this user guide, we will provide you with all the information you need to safely and effectively use your Violet or Blue ray device. We'll cover topics such as how to use the device, which electrode to use for different parts of the body, and what to expect during and after a treatment. We'll also provide you with tips for maintaining your device and troubleshooting any issues that may arise.

Before using your Violet or Blue ray device, it's important to read this guide thoroughly and follow all instructions carefully. While Violet or Blue ray therapy is generally considered safe, it's always a good idea to consult with a healthcare professional before using any new type of therapy. With proper use and care, your Violet or Blue ray device can be a valuable tool in promoting your overall health and wellness.

USER INSTRUCTIONS

Always read the user manual and follow all instructions carefully before using the device.

Clean the electrode before each use to ensure it is free of any debris or bacteria. Use a mild soap and water, or a disinfectant solution recommended by the manufacturer.

Choose the appropriate electrode for the area you wish to treat. The size and shape of the electrode can vary depending on the specific area being treated.

Apply a thin layer of conductive gel or cream to the skin before using the device. This will help ensure proper contact between the electrode and the skin, and may enhance the effectiveness of the treatment.

Start with the lowest intensity setting and gradually increase as needed. Do not exceed the recommended intensity level for the specific electrode being used.

Move the electrode in a slow, circular motion over the area being treated. Avoid staying in one spot for too long, as this can cause discomfort or skin irritation.

Treat each area for no more than 10-15 minutes at a time. Allow at least an hour between treatments to allow the body to rest and recover.

Drink plenty of water before and after each treatment to help flush toxins from the body.

Use caution when treating sensitive areas, such as the face. It's often best to start with a lower intensity setting and shorter treatment times in these areas.


Store the device in a cool, dry place and keep it out of reach of children. Follow any additional maintenance or storage instructions provided by the manufacturer.



Mushroom Wand

Most commonly used for any area of the body that requires relief from pain or inflammation.

Simply set to a comfortable frequency and use in circular motion while remaining in contact with skin. (Use water based gel to reduce friction on skin).



Y-Shape Wand

This wand can be used on neck area, arms or lower legs. Use up and down motion.


Can be used for pain or inflammation.



Comb Wand

The comb wand can be used on the scalp to stimulate hair growth. Use regular treatments for effect.

Use combing motion while keeping the comb teeth in contact with the scalp.



Tongue Wand

This wand can be used on the face to spot treat acne and pimples to reduce their appearance and severity quickly, by naturally killing the bacteria inside.



Rod Wand

This wand can be used for a general energy boost if you are feeling tired or lethargic.

Simply hold this wand in your hand and allow it to run for 5-10 minutes while your body gets a boost.



Bent Wand

This wand can again be used for acne, skin imperfections or small areas that require treatment.

Useful for getting into areas around nose and ears etc.

IMPORTANCE OF DIRECT SKIN CONTACT

Skin contact with a Violet or Blue Ray device is important because it allows the electrical energy from the device to be transmitted into the body, where it can have therapeutic effects. When the electrode of a Violet or Blue Ray device is in contact with the skin, the electrical current passes through the skin and into the underlying tissues. This can stimulate cellular activity, improve circulation, and promote healing.

In addition to allowing the electrical energy to be transmitted into the body, skin contact with a Violet or Blue Ray device is also important because it helps to prevent accidental electrical shocks. If the electrode of the device is not in contact with the skin, the electrical current may arc and cause a shock, which can be painful and potentially dangerous.

Another benefit of skin contact with a Violet or Blue Ray device is that it allows for targeted treatment of specific areas of the body. By placing the electrode directly on the skin, the electrical energy can be directed to a specific area, such as a painful joint or muscle, where it can have the greatest therapeutic effect.

It is important to note that skin contact with a Violet or Blue Ray device should be gentle and should not cause discomfort or pain. If the electrode is pressed too hard against the skin or moved too quickly, it can cause skin irritation or burns. It is also important to follow the manufacturer's instructions for use and to consult a healthcare professional if you have any concerns or questions about using a Violet or Blue Ray device.

BENEFITS OF USING A Violet or Blue RAY DEVICE

Pain relief:

Violet or Blue Ray therapy has been shown to help alleviate pain and discomfort caused by a variety of conditions, such as arthritis, muscle strains, and injuries.

Improved circulation:

The electrical current produced by the Violet or Blue Ray device may help improve blood circulation, which can lead to better oxygenation of tissues and organs.

Skin rejuvenation:

Violet or Blue Ray therapy is often used in cosmetic treatments to help reduce the appearance of fine lines, wrinkles, and acne. It can also help improve skin tone and texture.

Relaxation:

Violet or Blue Ray therapy may help promote relaxation by reducing muscle tension and promoting a sense of calm.

Respiratory relief:

The Violet or Blue Ray device can be used to treat respiratory conditions such as asthma, bronchitis, and sinusitis, by reducing inflammation and improving circulation in the affected areas.

Immune system support:

The electrical current produced by the Violet or Blue Ray device may help boost the immune system by stimulating white blood cell production and increasing overall immune function.

BENEFIT OF OZONE

Violet or Blue Rays emit a small amount of ozone gas when used. The ozone gas is produced when an electrical current is passed through the air between the electrode and the skin. Ozone has been shown to have several beneficial effects on the body, including:

Increased blood flow:

Ozone can help increase blood flow to the treated area by increasing the release of nitric oxide, a natural vasodilator. This can help improve circulation and promote healing.

Reduced inflammation:

Ozone has anti-inflammatory properties and can help reduce inflammation in the treated area. This can help reduce pain, swelling, and redness.

Improved immune function:

Ozone has been shown to stimulate the immune system and increase the production of white blood cells. This can help improve the body's ability to fight infections and promote healing.

Antimicrobial effects:

Ozone has antimicrobial properties and can help kill bacteria, viruses, and fungi. This can help prevent infections and promote healing.

COMMON PROTOCOLS

Skin Rejuvenation:

To promote skin rejuvenation, the Violet or Blue Ray electrode should be applied to the skin in a circular motion for 3-5 minutes per area, once or twice a week. This can help stimulate collagen production, improve skin tone and texture, and reduce the appearance of fine lines and wrinkles.

Pain Relief & Treating Inflammation:

To relieve pain, the Violet or Blue Ray electrode should be applied to the affected area for 3-5 minutes per session, 1-2 times a day. This can help reduce inflammation, increase blood flow, and promote tissue repair.

Hair Growth:

To promote hair growth, the Violet or Blue Ray electrode should be applied to the scalp in a combing or circular motion for 5-10 minutes per area, 1-2 times a week. This can help increase blood flow to the scalp, stimulate hair follicles, and promote hair growth.

CONTRAINDICATIONS

While Violet or Blue Ray therapy is generally considered safe for most people, there are some contraindications to be aware of. Here are some of the situations in which Violet or Blue ray therapy should not be used:

1. Pregnancy: Violet or Blue ray therapy should not be used during pregnancy, as it may affect the developing fetus.
2. Pacemaker: If you have a pacemaker or other implanted electrical device, you should not use a Violet or Blue ray device as it may interfere with the function of the device.
3. Metal implants: The electrical current produced by the Violet or Blue ray device may interfere with the function of metal implants, such as joint replacements or plates used to repair broken bones.
4. Cancer: If you have cancer or are undergoing cancer treatment, you should avoid using a Violet or Blue ray device without first consulting with a healthcare professional.
5. Epilepsy: Violet or Blue ray therapy may trigger seizures in people with epilepsy.
6. Skin conditions: If you have any open wounds, cuts, or infections on the skin, you should avoid using a Violet or Blue ray device as it may worsen the condition.
7. Sensitivity to electricity: Some people may be sensitive to electrical stimulation and may experience discomfort or pain during Violet or Blue ray therapy.
8. Medications: Certain medications, such as topical creams, may increase the risk of skin irritation or other side effects when used in conjunction with a Violet or Blue ray device. Check with a healthcare professional before using a Violet or Blue ray device if you are taking any medications.

As with any type of therapy, it's important to consult with a healthcare professional before using a Violet or Blue ray device, especially if you have any underlying health conditions or concerns.

DISCLAIMER

The use of a Violet or Blue ray device is not intended to diagnose, treat, cure, or prevent any medical condition. The information provided is for educational purposes only and is not a substitute for professional medical advice. Violet or Blue ray therapy should not be used as a replacement for any medical treatment prescribed by a healthcare professional. Always consult with a healthcare professional before using any new type of therapy, especially if you have any underlying health conditions or concerns. The use of a Violet or Blue ray device is at your own risk, and the manufacturer and distributor are not liable for any injury, damage, or adverse effects resulting from the use of this device. Keep the device out of reach of children and follow all instructions carefully for safe and effective use.