



Precision Vascular Access

Scripps

<https://www.testmenu.com/scripps>

Scripps Mission Valley Clinic 619-245-2840

LabCorp/VA

<https://www.labcorp.com/test-menu/search>

LabCorp Hillcrest 619-299-2831

VA Mission Valley Clinic 619-400-5032

CBC w/diff

Lavender top (EDTA)
1 ml whole blood



Lavender top (EDTA)

*Draw one tube for each test



CMP

*Gel green (lithium heparin plasma),
or Gold top (SST) acceptable

If using short gel green, submit one
tube for each test



Gold top (SST)



ESR

Lavender (EDTA) top
2 ml whole blood



Lavender top (EDTA)

2 ml whole blood

*Draw one tube for each test



CRP

*Gel green (lithium heparin plasma),
or Gold top (SST) acceptable








****MUST BE SPUN WITHIN 2 HOURS
OF COLLECTION****









If using short gel green, submit one
tube for each test



Gold top (SST)



<p>Trough, Vancomycin</p>	<p>Gold top (SST)</p> <p>**DRAW WITHIN 30 MINS OF NEXT DOSE**</p>		<p>Red top, plain (7 ml tube)</p> <p>*NO SST (gel)</p> <p>**DRAW WITHIN 30 MINS OF NEXT DOSE**</p>	
<p>Trough, Amikacin</p>	<p>Red top, plain (7 ml tube)</p> <p>*NO SST (gel)</p> <p>**DRAW WITHIN 30 MINS OF NEXT DOSE**</p>			
<p>Magnesium</p>	<p>*Gel green (lithium heparin plasma), or Gold top (SST) acceptable</p> <p>If using short gel green, submit one tube for each test</p>		<p>*Tiger Top (SST), or Gold top (SST) acceptable</p> <p>1 ml</p>	
<p>Phosphorus</p>	<p>*Gel green (lithium heparin plasma), or Gold top (SST) acceptable</p> <p>If using short gel green, submit one tube for each test</p>		<p>*Tiger Top (SST), or Gold top (SST) acceptable</p> <p>1 ml</p>	

<p>Triglycerides</p>	<p>*Gel green (lithium heparin plasma), or Gold top (SST) acceptable</p> <p>If using short gel green, submit one tube for each test</p>		<p>*Tiger Top (SST), or Gold top (SST) acceptable</p> <p>1 ml</p>	
<p>Prealbumin</p>	<p>*Gel green (lithium heparin plasma), or Gold top (SST) acceptable</p> <p>If using short gel green, submit one tube for each test</p>		<p>*Tiger Top (SST), or Gold top (SST) acceptable</p> <p>1 ml</p>	
<p>Lipase</p>	<p>Gold top (SST) 0.5 ml</p>		<p>Red top/Tiger top/Gold top or Green top (heparin) are all acceptable</p>	
<p>BNP</p>	<p>Gel green (lithium heparin plasma)</p>			
<p>Cryptococcal antigen</p>	<p>Gold top (SST) 1 ml</p>			

LFT/Hepatic
Function
Panel

Gel green (lithium heparin plasma)

****MUST BE SPUN WITHIN 2 HOURS
OF COLLECTION****



Red top/Tiger top/Gold top or
Green top (heparin) are all
acceptable

