

# Graf's Method cradles

The use of Graf's cradles is strongly recommended for scanning efficiently and delivering a correct diagnosis to patients, as mentioned in the reference manual **ESSENTIALS OF INFANT HIP SONOGRAPHY** (R. Graf, K. Lercher, S. Scott, T. Spieß) p. 25:

*"To perform quick and accurate scanning, a standard position for the infant and a standard scanning technique are needed. Used in combination with the cradle, the probe guide has proved very effective and facilitates scanning the infant hip, reduces scan time and considerably raising the level of precision.*

*By correctly positioning the infant and using the probe guide system, tilting errors can be virtually eliminated."*



**Legal disputes can occur if a child suffers long-term consequences due to delayed or inadequate treatment, as a result of diagnostic mistakes due to improper imaging techniques or incorrect interpretation of ultrasound images.**

## **Related products :**

### **Graf's Cradle (Rapid-Eco + Eco-Support)**

- Product code 5003
- Discount 10% for participants to courses
- size 27x51x80 cms - weight 7.3 KGS



### **Portable Graf's Cradle**

- Product code 5002
- Discount 10% for participants to courses
- size 44x38x22 cms - weight 4.8 KGS



**Contact customer service for enabling your discount:  
incoming@h-h-c.com or +41 22 8860003**