



LacTeck

# Pump2Baby Bottle



## WARNING

**Important: Please read all these instructions before use.**

- **Strangulation Hazard:** DO NOT allow tubing to wrap around baby's neck.
- ALWAYS use this product with adult supervision.
- Use bottle warmer or warm water from faucet to warm cold bottle. Do not microwave, as this may cause hotspots.
- Test temperature of warmed milk before giving to baby.
- Ensure baby is about half way between sitting up and lying down when feeding.
- Take product away from baby when baby is finished feeding.
- Do not use this product as a pacifier.
- This product is not a toy. Do not allow children to play with it.

## Cleaning:

- Prior to first use, sterilize all product components by boiling in water for five (5) minutes.
- After each use, wash each product component separately in warm soapy water. Or wash components separately on the top rack of dishwasher. Use a nipple brush (not included) to clean nipple.
- Use an infant syringe (not included) to clean the tube: fill the syringe with soapy warm water and push into the tubing. Use the syringe to rinse a couple times with clean water. Hang tube to dry.
- You can sterilize the parts by boiling water, microwave sterilizer, or steam sterilizer.

## Part Replacement:

Parts should be replaced for any of the following conditions:


- If tube becomes dirty or clogged, or if it no longer fits tightly on the posts.
- If nipple becomes cracked, sticky or damaged.
- If any parts are damaged or if leakage occurs.

- **CAUTION: To prevent a possible choking hazard**, test strength of nipple regularly by pulling hard on nipple. Inspect regularly, especially when baby is teething. If nipple becomes cracked, torn, or otherwise damaged, **discontinue use immediately**. Do not attempt to enlarge nipple hole, as this will cause tearing.
- **CAUTION: Pushing a nipple brush too firmly against the inside of the nipple can cause tearing of the nipple opening. Inspect the nipple opening regularly.**

## Assembly Instructions:

- Please ensure the tube is attached to the posts at both ends and be pushed all the way down onto each post to prevent leakage.
- The tube may be cut to any desired length to ensure a comfortable feeding position.
- Breast pump shield and storage caps used to convert to other bottles are not included. Please purchase separately from Medela and Dr. Brown's.

**Pump2Baby Bottle**



**Storage Bottle**



**Regular Bottle**



Please contact us at [contact@lacteck.com](mailto:contact@lacteck.com) for replacement parts or if you have any questions: