

Drug	Loading Dose	Maintenance Dose	Level
Midazolam	0.2 mg/kg IV q5min until seizure controlled or max 2 mg/kg	0.1-2.9 mg/kg/h infusion	N/A; titrate for EEG suppression
Propofol	2 mg/kg IV q5min until seizure controlled or max 10 mg/kg	2-15 mg/kg/h infusion	N/A; titrate for EEG suppression
Pentobarbital	5 mg/kg IV up to 50 min/min q5min until seizure controlled	1-10 mg/kg/h infusion	N/A; titrate for EEG suppression
Ketamine	1.5 mg/kg IV q5min until seizure controlled or max 4.5 mg/kg	1.2-7.5 mg/kg/h infusion	N/A; titrate for EEG suppression

DOWNLOAD: <https://tinurli.com/28u217>



7c23cee9bc

Related links:

- [hospital management project in java pdf download](#)
- [Torrent Download ArtCAM 2019 Portable](#)
- [Yistanita Duplicate Finder 3.9.6 Serial 93 mario uploaden shade](#)
- [Pirates Of The Caribbean On Stranger Tides Mp4 Download](#)
- [The 1 Was Framed full movie download in Italian](#)